

DECLARACIÓN DE CONFORMIDAD DE MATERIAL

EN CONTACTO CON ALIMENTOS

Alcance: Bandeja de pesada de aluminio.

Material: Lámina de aluminio con recubrimiento Softenol.

Referencia Labbox	Capacidad
WPAL-056-000	21 ml
WPAL-072-000	42 ml
WPAL-093-000	70 ml

Labbox Labware, bajo la responsabilidad del fabricante, declara que el producto anteriormente mencionado cumple con la legislación europea aplicable a los materiales y objetos destinados a entrar en contacto con alimentos, en particular:

- **Reglamento (CE) n.º 1935/2004** del Parlamento Europeo y del Consejo, de 27 de octubre de 2004, y sus modificaciones, sobre los materiales y objetos destinados a entrar en contacto con alimentos.
- **Reglamento (CE) n.º 2023/2006** sobre buenas prácticas de fabricación (BPF) para materiales y objetos destinados a entrar en contacto con alimentos.

Otras especificaciones:

- **Composición:** Lámina de aluminio con recubrimiento Softenol. Las posibles aleaciones utilizadas pertenecen a las series 8000 y 3000. El material puede contener un porcentaje de aluminio reciclado.
- **Alérgenos y OGM:** Los productos se declaran libres de alérgenos y de organismos genéticamente modificados (OGM). La documentación técnica hace referencia al Reglamento (UE) n.º 1169/2011 en relación con los alérgenos.
- **Condiciones de temperatura:** Rango de temperatura declarado de -40 °C a +400 °C.
- **Aplicaciones:** Aptos para uso en horno convencional, cocción al vapor, grill y congelación rápida. Compatibles con el sellado mediante film plástico.
- **Uso en microondas:** Aptos únicamente para determinados hornos microondas, bajo las condiciones de uso compatibles indicadas en el manual del equipo.
- **Resistencia química:** Buena resistencia al alcohol y a las grasas y aceites. Mala resistencia a ácidos, sustancias alcalinas y sales no diluidas.
- **Precauciones:** Evitar el contacto directo con llamas o resistencias y el contacto prolongado con alimentos fuertemente ácidos o salados.

Este documento ha sido generado electrónicamente y es válido sin firma: 22/10/2026

LABBOX GROUP

SPAIN • FRANCE • ITALY • GERMANY • NETHERLANDS

Main Office: Migjorn, 1 • 08338 Premià de Dalt (Barcelona) • Teléfono: + 34 937 077 970 • Fax: +34 937 909 532 • www.labbox.com

DECLARATION OF CONFORMITY FOR MATERIALS **IN CONTACT WITH FOOD**

Scope: Aluminium weighing dishes.

Material: Aluminium foil with Softenol coating.

Labbox Reference	Capacity
WPAL-056-000	21 ml
WPAL-072-000	42 ml
WPAL-093-000	70 ml

Labbox Labware, under the responsibility of the manufacturer, declares that the above-mentioned products comply with European legislation applicable to materials and articles intended to come into contact with food, in particular:

- **Regulation (EC) No 1935/2004** of the European Parliament and of the Council of 27 October 2004, and its amendments, on materials and articles intended to come into contact with food.
- **Regulation (EC) No 2023/2006** on Good Manufacturing Practice (GMP) for materials and articles intended to come into contact with food.

Other specifications:

- **Composition:** Aluminium foil with Softenol coating. Possible aluminium alloy grades belong to the 8000 and 3000 series. The material may contain a proportion of recycled aluminium.
- **Allergens and GMOs:** The products are declared free from allergens and genetically modified organisms (GMOs). The technical documentation refers to Regulation (EU) No 1169/2011 in relation to allergens.
- **Temperature conditions:** Declared temperature range from -40 °C to +400 °C.
- **Applications:** Suitable for conventional ovens, steam cooking, grilling and rapid freezing. Compatible with plastic film sealing.
- **Microwave use:** Suitable only for certain microwave ovens, under compatible operating conditions specified in the appliance instructions.
- **Chemical resistance:** Good resistance to alcohol, fats and oils. Poor resistance to acids, alkaline substances and undiluted salts.
- **Precautions:** Avoid direct contact with flames or heating elements and prolonged contact with strongly acidic or salty foods.

This document has been electronically generated and is valid without signature: 22/10/2026

LABBOX GROUP

SPAIN • FRANCE • ITALY • GERMANY • NETHERLANDS

Main Office: Migjorn, 1 • 08338 Premià de Dalt (Barcelona) • Teléfono: + 34 937 077 970 • Fax: +34 937 909 532 • www.labbox.com