

## GENERAL INFORMATION

**Product name :** Weighing pan, polystyrene, black

**Description :** Square or rhomboid shape

## TECHNICAL DATA

reference	capacity	h (mm)	l(mm)	w (mm)	type	pcs/pack
WPAN-044-250	7 ml	9	44	44	square shape	250
WPAN-080-250	100 ml	24	80	80	square shape	250
WPAN-135-250	250 ml	24	135	135	square shape	250
WPAN-058-250	5 ml	6	58	36	rhomboid shape	250
WPAN-059-250	30 ml	14	59	83	rhomboid shape	250
WPAN-099-125	100 ml	20	99	124	rhomboid shape	125

## PACKAGING AND LOGISTICS

reference	vol (l)	kg	TARIC	GTIN
WPAN-044-250	2,1	0,8	39269097	08434868083660
WPAN-080-250	6,8	2,5	39269097	08434868083677
WPAN-135-250	18,8	6	39269097	08434868083684
WPAN-058-250	1,6	0,5	39269097	08434868083691
WPAN-059-250	2,4	1,3	39269097	08434868083707
WPAN-099-125	3,2	2,9	39269097	08434868083714

## PRODUCT PHOTO



## INFORMACIÓN GENERAL

**Nombre del producto :** Bandeja de pesada de poliestireno, negra**Descripción :** Forma cuadrada o romboidal

## DATOS TÉCNICOS

referencia	capacidad	h (mm)	l (mm)	w (mm)	tipo	unidades por ref.
WPAN-044-250	7 ml	9	44	44	cuadrado	250
WPAN-080-250	100 ml	24	80	80	cuadrado	250
WPAN-135-250	250 ml	24	135	135	cuadrado	250
WPAN-058-250	5 ml	6	58	36	romboidal	250
WPAN-059-250	30 ml	14	59	83	romboidal	250
WPAN-099-125	100 ml	20	99	124	romboidal	125

## EMBALAJE Y DATOS LOGÍSTICOS

referencia	vol (l)	kg	TARIC	GTIN
WPAN-044-250	2,1	0,8	39269097	08434868083660
WPAN-080-250	6,8	2,5	39269097	08434868083677
WPAN-135-250	18,8	6	39269097	08434868083684
WPAN-058-250	1,6	0,5	39269097	08434868083691
WPAN-059-250	2,4	1,3	39269097	08434868083707
WPAN-099-125	3,2	2,9	39269097	08434868083714

## FOTO DEL PRODUCTO



## INFORMATIONS GÉNÉRALES

**Nom produit :** Nacelle de pesée en polystyrène, noire**Description :** Forme carrée ou en losange

## DONNÉES TECHNIQUES

référence	capacité	h (mm)	long. (mm)	larg. (mm)	type	unités par ref.
WPAN-044-250	7 ml	9	44	44	carré	250
WPAN-080-250	100 ml	24	80	80	carré	250
WPAN-135-250	250 ml	24	135	135	carré	250
WPAN-058-250	5 ml	6	58	36	losange	250
WPAN-059-250	30 ml	14	59	83	losange	250
WPAN-099-125	100 ml	20	99	124	losange	125

## EMBALLAGE ET LOGISTIQUE

référence	vol (l)	kg	TARIC	GTIN
WPAN-044-250	2,1	0,8	39269097	08434868083660
WPAN-080-250	6,8	2,5	39269097	08434868083677
WPAN-135-250	18,8	6	39269097	08434868083684
WPAN-058-250	1,6	0,5	39269097	08434868083691
WPAN-059-250	2,4	1,3	39269097	08434868083707
WPAN-099-125	3,2	2,9	39269097	08434868083714

## PHOTO PRODUIT



## INFORMAZIONE GENERALE

**Nome del prodotto** : Navicella di pesata in polistirene, nero**Descrizione** : Forma quadrata o romboidale

## DATI TECNICI

referenza	capacità	h (mm)	L (mm)	w (mm)	tipo	unità per ref.
WPAN-044-250	7 ml	9	44	44	quadrato	250
WPAN-080-250	100 ml	24	80	80	quadrato	250
WPAN-135-250	250 ml	24	135	135	quadrato	250
WPAN-058-250	5 ml	6	58	36	romboidale	250
WPAN-059-250	30 ml	14	59	83	romboidale	250
WPAN-099-125	100 ml	20	99	124	romboidale	125

## IMBALLAGGIO E DATI LOGISTICI

referenza	vol (l)	kg	TARIC	GTIN
WPAN-044-250	2,1	0,8	39269097	08434868083660
WPAN-080-250	6,8	2,5	39269097	08434868083677
WPAN-135-250	18,8	6	39269097	08434868083684
WPAN-058-250	1,6	0,5	39269097	08434868083691
WPAN-059-250	2,4	1,3	39269097	08434868083707
WPAN-099-125	3,2	2,9	39269097	08434868083714

## FOTO DEL PRODOTTO

