

PCR 8-TUBE STRIPS

PCRS-0A8-120

astik'S


Date: 31-03-2022

Material: Polypropylene USP VI

Colour: Natural, transparent

Capacity: 0,2 ml

Features: Caps attached to each tube. Human DNA free and certified DNase, RNase, and PCR inhibitors free



PCR 8-TUBE STRIPS

PCRS-0S8-120

astik'S

Date: 31-03-2022

Material: Polypropylene USP VI

Colour: Natural, ransparent

Capacity: 0,2 ml

Features: Caps included, separately. Human DNA free and certified DNase, RNase, and PCR inhibitors free