

qPCR 8-TUBE STRIPS

PCRS-8C1-120

astik'S


Date: 31-03-2022

Material: Polypropylene USP VI

Colour: Natural, frosted

Capacity: 0,1 ml

Features: Optical caps included, separately. Human DNA free and certified DNase, RNase, and PCR inhibitors free



qPCR 8-TUBE STRIPS

PCRS-8C2-120

astik'S


Date: 31-03-2022

Material: Polypropylene USP VI

Colour: Natural, transparent

Capacity: 0,2 ml

Features: Optical caps included, separately. Human DNA free and certified DNase, RNase, and PCR inhibitors free



qPCR 8-TUBE STRIPS

PCRS-8W1-120

astik'S


Date: 31-03-2022

Material: Polypropylene USP VI

Colour: White, opaque

Capacity: 0,1 ml

Features: Optical caps included, separately. Human DNA free and certified DNase, RNase, and PCR inhibitors free



qPCR 8-TUBE STRIPS

PCRS-8W2-120

astik'S


Date: 31-03-2022

Material: Polypropylene USP VI

Colour: White, opaque

Capacity: 0,2 ml

Features: Optical caps included, separately. Human DNA free and certified DNase, RNase, and PCR inhibitors free



qPCR 8-TUBE STRIPS

PCRS-8U1-120

astik'S

Date: 08-03-2023

Material: Polypropylene USP VI

Colour: Ultra clear, transparent

Capacity: 0,1 ml

Features: Optical caps included, separately. Human DNA free and certified DNase, RNase, and PCR inhibitors free